

3DVTM01-3.0ML



The 3D MedSupply and Manufacturing VTM Tube is made in our ISO 5 clean rooms here in the USA with FDA Registration # 3020027409. We follow the strict CDC and FDA guidelines during the manufacturing process. The VTM is designed to facilitate the identification of all Viruses including, SARS-CoV-2, Influenza A/B, RSV, as well as chlamydia and mycoplasma (bacteria) with any PCR technology by providing:

- Safe transport and preservation of the specimen
- Verified nucleic acid extraction and qRT-PCR, RT-PCR, qPCR, and RT-qPCR



PRINCIPLES

The 3DVTM01-3.0ML Transit Tube contains a viral transport medium (VTM; Culture Media, Non-Propagating, Transport) intended to be inoculated with a Nasopharyngeal (NP), Oropharyngeal (OP), or Anterior Nares (AN) synthetic fiber swab specimens. It is then transported appropriately to the lab and analyzed with validated qRT-PCR, RT-PCR, qPCR and RT-qPCR assays for the detection of all Viruses, including severe acute respiratory syndrome (SARS-Cov-2), Influenza A/B, RSV, as well as chlamydia and mycoplasma on humans.

SPECIFICITY

3DVTM01-3.0ML suppresses the growth of other bacteria and fungi that may be present in clinical samples from the human respiratory system.

SHELF LIFE

Can be stored at 2-8°C for up to 12 months and room temperature for 4 months without deterioration of performance, but not past the expiration date on the tube.

SPECIMEN

Nasopharyngeal, (NP), oropharyngeal (OP), or Anterior Nares (AN) swab

AVAILABILITY

Contact us for availability.

TRANSPORT

Inoculated tubes should be transported within 72 hours after inoculation and maintained at 2-8°C



Establishment Registration & Device Listing

| Establishment Name | Registration Number | Current Registration Yr |
|--|---------------------|-------------------------|
| 3D MedSupply & Manufacturing, LLC | FL/USA 3020027409 | 2022 |
| Culture Media, Non-Selective And Non-Differential - 3DSWB01 | | Manufacturer |
| Transport Medium, Notified Per The Vtm Guidance - 3DVTM01-1.5ML, 3DVTM01-3.0ML | | Manufacturer |



IN STOCK



Box Contents

- (50) Sterile 3ml VTM tubes
- (50) Sterile Swabs
- (50) Biohazard Bags



TO ORDER NOW CALL

Office (866) 30-3D-MED / (866) 303-3633
 Fax (855) 96-3D-MED / (855) 963-3633
 Or email: Info@3dmedsupply.com

4303 Vineland Rd. Suite F8
 Orlando, FL 32811
 www.3DMedicalSupply.com

